

# Chloride Ion

ELIT 8261 · ELIT Ion Selective Electrode · Anion

Cl<sup>-</sup>

>> View full ion guide: [elyxir.co.uk/resources/ions/chloride](http://elyxir.co.uk/resources/ions/chloride)

Order this electrode: [elyxir.co.uk/products](http://elyxir.co.uk/products)

## PHYSICAL SPECIFICATIONS

Body Length	130 mm (excl. contact) / 140 mm (incl.)
Body Diameter	8 mm
DC Resistance (25°C)	< 0.5 MOhm
Min. Sample Volume	5 ml

## ELECTRODE SPECIFICATIONS

Electrode Model	ELIT 8261
Ion	Chloride (Cl <sup>-</sup> )
Ion Type	Anion
Valence	1
Membrane Type	Solid-state poly-crystalline membrane
Molar Mass	35.453 g/mol
1000 ppm equiv.	0.028 M

## OPERATIONAL PARAMETERS

Preconditioning	1000 ppm Chloride standard
Preconditioning Time	Min. 5 minutes
Detection Range	1 to 35,000 ppm ( $3 \times 10^{-5}$ to 1 M)
Electrode Slope	$54 \pm 5$ mV/decade
pH Range	pH 1 to 12
Temperature Range	0 to 80 °C
Response Time	< 10 seconds (90% response)
Potential Drift	< 3 mV/day in 1000 ppm (8 hours)

## REAGENTS & STANDARDS

Reference Electrode	Double junction (ELIT 003). Outer filling solution: 0.1M CH <sub>3</sub> COOLi.
ISAB / Buffer	5M NaNO <sub>3</sub> — Add 2% v/v.
Standard Prep	Dissolve 1.649 g anhydrous sodium chloride (NaCl) in 1 litre deionised water.

## TYPICAL APPLICATIONS

- Food & Beverage Analysis
- Water Quality Monitoring
- Agriculture & Soil Analysis
- Industrial Process Control
- Environmental Monitoring

## SELECTIVITY COEFFICIENTS (INTERFERENCE DATA)

Interfering Ion	Selectivity Coeff.	Note
Iodide (I <sup>-</sup> )	>1	CRITICAL: the Cl membrane is far more sensitive to I <sup>-</sup> than to Cl <sup>-</sup> . High iodide will irreversibly damage the membrane.
Bromide (Br <sup>-</sup> )	high	Strong interference — must be absent or at insignificant concentrations.
Cyanide (CN <sup>-</sup> )	high	Strong interference. Do NOT add acid to cyanide-containing samples — risk of lethal HCN gas.
Sulphide (S <sup>2-</sup> )	high	All poly-crystalline (AgS-based) membranes are unreliable in presence of Ag or S ions.
Silver (Ag <sup>+</sup> )	high	—

SC = approximate apparent increase in measured concentration caused by 1 unit of interferent. Error% = ((interferent conc × SC) / target conc) × 100.

## CALIBRATION & SAMPLE PREPARATION

---

Calibrate with 1000, 100, 10, 1 ppm Cl solutions. For high ionic strength samples (> 0.01 M), add 2 ml ISAB to each 100 ml standard.

Low ionic strength: immerse in 50–100 ml sample. High ionic strength: add 2 ml ISAB to 100 ml sample and stir well.

## SPECIAL ANALYTICAL PROCEDURES

---

### Butter

Heat ~10 g butter with 100 ml dilute HNO<sub>3</sub> until fat melts. Stir to extract chloride. Cool until fat solidifies. Pipette 25 ml of solution through fat layer and dilute 1:1 with deionised water. Standards must be mixed 1:1 with dilute HNO<sub>3</sub> ISAB.

### Milk

Mix 1 volume milk with 2 volumes dilute HNO<sub>3</sub> (60 ml conc. HNO<sub>3</sub> per litre) and shake for 1 minute. Nitric acid acts as ISAB. Measure by Direct Potentiometry. Add HNO<sub>3</sub> ISAB to all standards in 2:1 ratio.

### Meat

Liquidise ~1 g meat with 100 ml water. Filter, mix 25 ml with 25 ml of 0.5M HNO<sub>3</sub> ISAB. Result: ppm × 100 ÷ sample weight = µg/g.

## ANALYTICAL NOTES

---

- The chloride electrode will only give reliable results if iodide, bromide, cyanide, sulphide, and silver are absent or at insignificant levels.
- A sodium bromate buffer can remove up to 1000 ppm Br or I and 500 ppm sulphide. Prepare and use in a well-ventilated area.

## SAFETY & HAZARDS

---

**! Do NOT add sodium bromate buffer to samples containing significant cyanide — the acid solution may liberate lethal HCN gas.**  
**! Irreversible membrane damage will result from exposure to high iodide concentrations.**

---

This document is provided for guidance only. Specifications subject to change without notice. For technical support contact sales@nico2000.net or call 020 8422 6779.