

# Perchlorate Ion

ELIT 8061 · ELIT Ion Selective Electrode · Anion

ClO<sub>4</sub><sup>-</sup>

>> View full ion guide: [elyxir.co.uk/resources/ions/perchlorate](https://elyxir.co.uk/resources/ions/perchlorate)

Order this electrode: [elyxir.co.uk/products](https://elyxir.co.uk/products)

## PHYSICAL SPECIFICATIONS

Body Length	130 mm (excl. contact) / 140 mm (incl.)
Body Diameter	8 mm
DC Resistance (25°C)	< 2.5 MOhm
Min. Sample Volume	5 ml

## ELECTRODE SPECIFICATIONS

Electrode Model	ELIT 8061
Ion	Perchlorate (ClO <sub>4</sub> <sup>-</sup> )
Ion Type	Anion
Valence	1
Membrane Type	Solid-state PVC polymer matrix membrane
Molar Mass	99.451 g/mol
1000 ppm equiv.	0.010 M

## OPERATIONAL PARAMETERS

Preconditioning	1000 ppm Perchlorate standard
Preconditioning Time	Min. 5 minutes
Detection Range	0.2 to 9,960 ppm (2×10 <sup>-6</sup> to 0.1 M)
Electrode Slope	54 ± 5 mV/decade
pH Range	pH 0 to 11
Temperature Range	0 to 50 °C
Response Time	< 10 seconds (90% response)
Potential Drift	< 3 mV/day in 1000 ppm (8 hours)

## REAGENTS & STANDARDS

Reference Electrode	Single junction (ELIT 001).
ISAB / Buffer	1M CH <sub>3</sub> COONa (sodium acetate) — Add 2% v/v.
Standard Prep	Dissolve 1.231 g anhydrous sodium perchlorate (NaClO <sub>4</sub> ) in 1 litre deionised water.

## TYPICAL APPLICATIONS

- Water Quality Monitoring
- Environmental Monitoring
- Industrial Process Control
- Research

## SELECTIVITY COEFFICIENTS (INTERFERENCE DATA)

Interfering Ion	Selectivity Coeff.	Note
Thiocyanate (SCN <sup>-</sup> )	0.03	—
Iodide (I <sup>-</sup> )	0.02	Tolerated up to the same concentration as perchlorate (~2% enhancement).
Nitrate (NO <sub>3</sub> <sup>-</sup> )	0.02	Tolerated up to the same concentration as perchlorate (~2% enhancement).
Chloride (Cl <sup>-</sup> )	0.0003	Can be present up to 100x more concentrated before significant error.
Phosphate (PO <sub>4</sub> <sup>3-</sup> )	0.0002	—
Acetate (CH <sub>3</sub> COO <sup>-</sup> )	0.0001	—

SC = approximate apparent increase in measured concentration caused by 1 unit of interferent. Error% = ((interferent conc × SC) / target conc) × 100.

## CALIBRATION & SAMPLE PREPARATION

---

Calibrate with 1000, 100, 10, 1, 0.1 ppm  $\text{ClO}_4$  solutions. For high ionic strength samples ( $> 0.01 \text{ M}$ ), add 2 ml ISAB to each 100 ml standard.

Low ionic strength: immerse in 50–100 ml sample. High ionic strength: add 2 ml ISAB to 100 ml sample and stir well.

## ANALYTICAL NOTES

---

- Wide pH operating range (0 to 11).
- Nitrate can be tolerated at equal concentration to perchlorate (causes ~2% enhancement).
- Chloride can be present up to 100x the perchlorate concentration before causing significant error.

---

This document is provided for guidance only. Specifications subject to change without notice. For technical support contact [sales@nico2000.net](mailto:sales@nico2000.net) or call 020 8422 6779.